

Sentina® Ivory Latex

Powder-free examination glove



Sentina® Ivory Latex Powder-free examination glove



Product information:

- High wearing comfort and excellent tactile sensitivity thanks to optimised grammage.
- Secure grip thanks to textured fingertip.
- With synthetic inner coating for quick and easy donning.

Fields of application:

- For use in hospitals and in outpatient settings.
- For professional cleansing.
- Suitable for laboratories, pharmaceutical companies and for contact with food.

Sentina® Ivory Latex

Powder-free examination glove,
non-sterile, off-white

Size	REF	Shipping units (pieces)
XS	163 040	100/1,000
S	163 041	100/1,000
M	163 042	100/1,000
L	163 043	100/1,000
XL	163 044	90/900

Product features:

Product type:	Powder-free examination and disposable protective gloves, non-sterile
Material:	Natural rubber latex
Colour:	Off-white
Glove shape:	Flat with rolled edge, ambidextrous
Glove surface:	Textured fingertip
Size/overall length:	XS, S, M, L, XL / median 240 mm (according to EN 455-2)
Wall thickness: (Simple wall thickness)	Palm: min. 0.08 mm, typ. 0.10 mm Finger: min. 0.08 mm, typ. 0.12 mm
Barrier properties:	AQL 1.5 (according to EN 455-1)
Tear resistance:	Median ≥ 6 newtons (according to EN 455-2)
Shelf life:	5 years (in original packaging when stored in accordance with ISO 2230)
Packaging:	100/90 pcs per dispenser box, 10 dispenser boxes per transport carton
Labelling and conformity to standards:	EN 21420, EN ISO 374-1, EN ISO 374-2, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ASTM F1671, EN 16523-1
Intended purpose:	Disposable medical gloves class I according to Medical Device Regulation (EU) 2017 / 745. Disposable protective gloves category III (temporary protection against chemical effects)* according to PPE regulation (EU) 2016 / 425. Suitable for contact with food according to regulation (EC) 1935 / 2004.

Additional information: This product contains natural rubber latex and chemical accelerators.

In case of intolerance, do not use these gloves.

A list of the substances contained in the glove is available on request.

* Please note that the product must be used in accordance with its chemical resistance and the respective declaration of conformity. Failure to observe this information, in particular with regard to chemical resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Lohmann & Rauscher is not liable for incorrect use of the gloves nor will Lohmann & Rauscher be liable for improper storage and handling. In case of indications of damage or unusual signs of wear and tear to the gloves, the gloves should not be used. In case of any doubt, expert advice should be sought before use. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments.

Publisher

Lohmann & Rauscher GmbH & Co. KG, Irlicher Straße 55, D-56567 Neuwied
Lohmann & Rauscher GmbH, Kirchengasse 17, A-2525 Schönaun an der Triesting

Contact address

Lohmann & Rauscher GmbH & Co. KG, Westerwaldstraße 4, D-56579 Rengsdorf